

***Technical Publication***



**DICOM CONFORMANCE STATEMENT**

**BrainLAB DICOM Service 3.50  
(Intel)**

**Revision 1.0**

## Table of contents

<b>INTRODUCTION</b> .....	<b>3</b>
PURPOSE.....	3
ABBREVIATIONS.....	3
REFERENCES.....	3
<b>IMPLEMENTATION MODEL</b> .....	<b>4</b>
APPLICATION DATA FLOW DIAGRAM.....	4
FUNCTIONAL DEFINITION OF APPLICATION ENTITY (AE).....	5
SEQUENCING OF REAL WORLD ACTIVITIES.....	5
<b>APPLICATION ENTITY SPECIFICATIONS</b> .....	<b>6</b>
STORAGE AE SPECIFICATIONS.....	6
ASSOCIATION ESTABLISHMENT POLICIES.....	7
<i>General</i> .....	7
<i>Number of associations</i> .....	7
<i>Asynchronous nature</i> .....	7
<i>Implementation identifying information</i> .....	7
<i>Association acceptance by real-world activity</i> .....	7
<i>Real-world activity for Receive message operations</i> .....	8
<i>Associated real-world activity for Receive message operations</i> .....	8
<i>Presentation context acceptance criterion for Receive message operations</i> .....	8
<i>Transfer syntax selection policies for Receive message operations</i> .....	8
<b>COMMUNICATION PROFILES</b> .....	<b>9</b>
SUPPORTED COMMUNICATION STACKS.....	9
TCP/IP STACK.....	9
PHYSICAL MEDIA SUPPORT.....	9
POINT TO POINT STACK.....	9
<b>EXTENSIONS/ SPECIALIZATIONS/ PRIVATIZATIONS</b> .....	<b>10</b>
STANDARD EXTENDED/ SPECIALIZED/ PRIVATE SOP'S.....	10
PRIVATE TRANSFER SYNTAXES.....	10
<b>CONFIGURATION</b> .....	<b>10</b>
AE TITLE/ PRESENTATION ADDRESS MAPPING.....	10
CONFIGURABLE PARAMETERS.....	10
<b>SUPPORT OF EXTENDED CHARACTER SETS</b> .....	<b>10</b>

## ***Introduction***

### ***Purpose***

This is a conformance statement for the DICOM extension of the BrainLAB product. The extension (Windows NT Service) supports DICOM storage services as a Service Class Provider (SCP).

This DICOM Conformance Statement is written according to part PS 3.2 of [1].

### ***Abbreviations***

ACR	American College of Radiation
AE	DICOM Application Entity
NEMA	National Electrical Manufactures Association
SCU	DICOM Service Class User
SCP	DICOM Service Class Provider
SOP	Service Object Pair
UID	Unique Identifier
VR	Value Representation

### ***References***

[1] Digital Imaging and Communications in Medicine (DICOM) 3.0,  
NEMA PS 3.1-3.13



The BrainLAB DICOM implementation is based on the MergeCOM-3 DICOM library from Merge Healthcare.

***Implementation Model***

The BrainLAB DICOM Application Entity is an implementation of a DICOM Storage Service Class Provider (SCP) which can receive DICOM messages from a DICOM Storage Service Class User (SCU).

The AE Title of the BrainLAB storage SCP is configurable by a configuration file. But for the reason of simplification, the AE Title will be named BRAINLAB\_SCP further in the text.

***Application Data Flow Diagram***




Figure 2A: BRAINLAB\_SCP application data flow diagram

***Functional definition of Application Entity (AE)***

All communications and message transfer with the remote application is accomplished utilizing the DICOM protocol over a network using the TCP/ IP protocol stack.

The BRAINLAB\_SCP will respond, if asked, with the Verification SOP Class UID as a SCP for one of its implemented SOP classes.

BRAINLAB\_SCP waits for an association to accept at the TCP/ IP port number that is configured at the time this application is initiated. When an association request is received with valid connection criteria, BrainLAB SCP responds with a list of SOP class UIDs that it will accept. It then waits for a Store request. If a Store is received, then all incoming messages that are conformant to the association are written to files on disk.

***Sequencing of real World Activities***

Not applicable.

## Application Entity Specifications

### Storage AE specifications

BRAINLAB\_SCP provides Standard Conformance to the following DICOM v3.0 SOP Class as a Verification Service Class Provider (SCP). As a SCP it sends out an Echo response after it receives an Echo request from a remote AE.

SOP Class UID	SOP Class Name
1.2.840.10008.1.1	Verification SOP Class

Table 3A: Valid SCP Verification SOP Class for the BrainLAB SCP AE

BRAINLAB\_SCP provides Standard Conformance to the following DICOM v3.0 Service Object Pair (SOP) Classes as a Storage Service Class Provider (SCP).

SOP Class UID	SOP Class Name
MR Image Storage	1.2.840.10008.5.1.4.1.1.4
CR Image Storage	1.2.840.10008.5.1.4.1.1.1
CT Image Storage	1.2.840.10008.5.1.4.1.1.2
Stand-alone Curve Storage	1.2.840.10008.5.1.4.1.1.9
Stand-alone Modality LUT Storage	1.2.840.10008.5.1.4.1.1.10
Stand-alone Overlay Storage	1.2.840.10008.5.1.4.1.1.8
Secondary Capture (SC) Image Storage	1.2.840.10008.5.1.4.1.1.7
Ultrasound (US) Image Storage	1.2.840.10008.5.1.4.1.1.6.1
US Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1
Stand-alone VOI LUT Storage	1.2.840.10008.5.1.4.1.1.11
Nuclear Medicine (NM) Image Storage	1.2.840.10008.5.1.4.1.1.20
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1
X-Ray Angiographic Bi-plane Image Storage	1.2.840.10008.5.1.4.1.1.12.3
X-Ray Radiofluoroscopic (RF) Image Storage	1.2.840.10008.5.1.4.1.1.12.2
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3
RT Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.4
RT Brachy Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.6
RT Treatment Summary Record Storage	1.2.840.10008.5.1.4.1.1.481.7
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128
Standalone PET Curve Storage	1.2.840.10008.5.1.4.1.1.129
Digital X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.1.1
Digital X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1.1
Digital Intra-oral X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.3
Digital Intra-oral X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.3.1
Digital Mammography Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.2
Digital Mammography Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2
VL Photographic Image	1.2.840.10008.5.1.4.1.1.77.1.4
VL Slide-Coordinates Microscopic Image	1.2.840.10008.5.1.4.1.1.77.1.3
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6

Table 3B: Valid SCP Storage SOP Classes for the BrainLAB SCP AE

### ***Association establishment policies***

#### ***General***

The BRAINLAB\_SCP application will wait for an association as a Storage SCP. When a Store request is received, the corresponding messages are saved to files on disk.

The maximum PDU size is 28672 bytes.

#### ***Number of associations***

The BRAINLAB\_SCP AE allows multiple simultaneous Store associations. The maximum number of simultaneous associations is 5.

#### ***Asynchronous nature***

BRAINLAB\_SCP does not support asynchronous communication (multiple outstanding transactions over a single association).

#### ***Implementation identifying information***

The Implementation Class Unique Identifier (UID) for the BRAINLAB\_SCP Application Entity is: 1.2.276.0.20.1.1

#### ***Association acceptance by real-world activity***

The BRAINLAB\_SCP client application accepts an association for the appropriate Storage Service Class that corresponds to the set of messages requested to be transferred. The association is closed by the Storage Service Class user, which initiated the association.

BRAINLAB\_SCP is able to abort the association when an error occurs.

**Real-world activity for Receive message operations**

BRAINLAB\_SCP waits for an association and offers to do the Storage Service. The association is closed after an error or when the initiator requests that it be closed.

**Associated real-world activity for Receive message operations**

Once the association has been established, the BRAINLAB\_SCP waits for transmission of conformant Storage Service messages.

Presentation Context Table				
Abstract Syntax	Transfer Syntax		Role	Extended Negotiation
<b>All Storage SOP Classes declared above in Table 3B</b>	Name List	UID List	SCP	None
	DICOM Implicit VR Little Endian	1.2.840.10008.1.2		
	DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2		
	DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1		

Table 3C: Receive Message Presentation Contexts of BRAINLAB\_SCP

**SOP specific conformance for all storage SOP Classes**

No known SOP specific conformance issues.

**Presentation context acceptance criterion for Receive message operations**

BRAINLAB\_SCP will accept the verification or storage SOP classes that are listed above. In the event of receiving an unknown presentation context the BRAINLAB\_SCP will reject the association request.

**Transfer syntax selection policies for Receive message operations**

BRAINLAB\_SCP supports the Implicit VR Little Endian, the Explicit VR Little Endian and the Explicit VR Big Endian transfer syntax. Any proposed presentation context, which includes one of these transfer syntaxes, will be accepted. Any proposed presentation context that does not include one of these transfer syntaxes will be rejected.



## ***Communication profiles***

### ***Supported Communication Stacks***

BRAINLAB\_SCP supports the DICOM upper layer using TCP/IP.

#### ***TCP/IP Stack***

The TCP/IP stack is inherited from the Windows® NT™ Operating system.

#### ***Physical Media Support***

Ethernet v2.0, IEEE 802.3

#### ***Point to Point Stack***

The 50-pin ACR-NEMA connection is not applicable to this product.

### ***Extensions/ specializations/ privatizations***

#### ***Standard extended/ specialized/ private SOP's***

None supported.

#### ***Private Transfer Syntaxes***

None supported.

### ***Configuration***

The BRAINLAB\_SCP application references one configuration file, which is located in the operating systems directory.

#### ***AE title/ presentation address mapping***

Presentation address mapping is configured in the BrainLAB initialization file. The presentation address of a SCP application as a provider is specified by configuring the listen Port in the BrainLAB initialization file, and specifying the AE title.

#### ***Configurable parameters***

The BrainLAB initialization file can be used to set the basic configuration parameters. This includes the TCP/IP listen port, time-out and the AE Title.

### ***Support of extended character sets***

Not supported.

Last page of document

BrainLAB AG 2001